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CONTENTS

1	General.....	27
1.1	EBS Article numbers	27
1.2	Transport damage / Unpacking the Device	27
1.3	Scope of guarantee	28
1.4	Contact Address / Manufacturer	28
2	Safety	29
2.1	Intended use.....	29
2.2	Microbiological examination of decentralized disinfectant dispensers	29
2.3	General Safety Information	30
2.4	Maintenance and Repair	30
2.5	Safety measures (to be performed by the owner/operator)	30
2.5.1	Trained person	31
2.5.2	Qualified person	31
2.6	Responsibilities of the Owner.....	31
2.7	System Separation	31
2.8	Emphases / Numberings	32
2.9	Safety symbols	32
2.10	Machine identification	32
3	Design.....	33
4	Functional Description and Operation.....	34
4.1	Description of Display Screen and Operating Controls.....	34
5	Connection and Installation	35
5.1	Wall Mounting.....	35
5.2	Water connection	36
5.3	Electrical Connection.....	36
5.4	Shifting Suction Pipe and Mains Cable to Left-Hand Side of Housing	37
6	Settings and Start-Up.....	38
6.1	Start-Up	38
6.2	Setting the Concentration.....	39
6.2.1	Checking the Concentration Setting.....	40
6.2.2	Checking the Safety Equipment	40
6.2.3	Functional Check of the Level Suction Pipe	40
6.2.4	Functional Check of the Flow Control Device.....	41
6.2.5	Changing over to a New Product.....	41
7	Troubleshooting.....	42
8	Maintenance and Repair Information.....	43
8.1	Microbiological examination of decentralized disinfectant dispensers.....	43
8.2	General maintenance work	43
8.2.1	Maintenance Table.....	44
8.3	Recommendations for the replacement of wearing parts	44
8.4	Repair Instructions	44
9	Technical Data.....	45
9.1	General Data	45
9.2	Connection Diagram.....	46
9.3	Connection Material and Accessories.....	47
10	Dismantling / Disposal.....	48
11	Declaration of Conformity	49

1 General

This Operating Manual contains all the information required for the installation, start-up, and maintenance of the **DG 1-II** disinfectant metering unit.

The relevant chapters of this manual must be read and noted before installing the device, using it for the first time, and before carrying out any maintenance or repair work.

The **DG 1-II** is a metering unit for adding disinfectant in order to prepare a ready-to-use disinfecting solution.

DG1- II / EPDM, Article No. 102410:

Metering pump, Flow control and Suction pipe with EPDM O-Rings.
Recommended use in connection with products **without glucoprotamine**.

DG1- II / SILICONE, Article No. 1024:

Metering pump, Flow control and Suction pipe with SILOCONE O-Rings.
Recommended use in connection with products **containing glucoprotamine**.

NOTE

This Operating Manual is a component of the equipment and must be available at all times to the operating and maintenance staff. The safety information contained therein must be observed. In the event the equipment is resold, the Operating Manual must always be supplied with it.

The German chapters of this manual constitute the **ORIGINAL OPERATING INSTRUCTIONS** and take legal precedence. All other languages are translations of the **ORIGINAL OPERATING INSTRUCTIONS**.

If you have questions, please contact us as outlined in chapter **1.4** "Contact address / manufacturer".

1.1 EBS Article numbers

Within this operating instruction article numbers, and EBS numbers are represented. EBS numbers are **Ecolab** internal article numbers and used only "concern internal".

1.2 Transport damage / Unpacking the Device

When removing the **DG 1-II** from its box, make sure that you do not rest the device on its lower side, since this can damage the dosing pipe that protrudes from the connection for the swivel spout.

CAUTION

If transport damage is discovered at the time of unpacking, the **DG1-II** unit may not be installed! The carrier in question must be notified immediately and a damage

ATTENTION

confirmation slip must be requested. A copy of the damage confirmation slip must be sent to **Ecolab Engineering GmbH** immediately.

1.3 Scope of guarantee

The manufacturer only guarantees the operational safety, reliability, and performance of the **DG1-II** under the following conditions:

- The installation, connection, setting, maintenance, and repair work is carried out by qualified and authorised specialists.
- The metering unit must be used in accordance with the specifications as given in the technical manual.
- Only original equipment spare parts are used for repairs.
- Only “validated” Ecolab-disinfectants are used. If other than Ecolab products are used, it is recommended that the disinfectant is validated for the dosing unit, according to BAM 2004-47;67-72.

For the remainder the general warranty and service terms of **Ecolab Engineering GmbH** shall apply.

1.4 Contact Address / Manufacturer

Ecolab Engineering GmbH

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D-83313 Siegsdorf

Telephone (+49) 86 62 / 61 0

Fax (+49) 86 62 / 61 2 35

eMail: engineering-mailbox@ecolab.com

2 Safety

**ATTENTION**

The **DG1-II** system must be used with validated products
No warranty for the unit function when using unvalidated products.

**CAUTION**

Safety warnings and emphases must be observed in all cases!

**ATTENTION**

2.1 Intended use

The **DG 1-II** is a mixing equipment for disinfectants.

It is used for the production of a ready-to-use disinfectant solution and has been developed, designed and built exclusively for the industrial and commercial use. A private use of the equipment is excluded.

It is part of the intended use of dosing system to follow the manufacturer's instructions for maintenance, service, repair and disposal.

The service life of the equipment according to the properly conducted maintenance work is approximately 10 years . An audit (potentially followed by a general overhaul) is required subsequently by the manufacturer or a specialist company.

2.2 Microbiological examination of decentralized disinfectant dispensers

The end user is responsible for the proper operation of the dosing device according to the recommended process from the RKI "requirements for design, features and operation of decentralized disinfectant dispensers "(Federal Health Bulletin 2004 47: 67-72.).

In the device manual the annual required hygiene-technical checks are described in section related to maintenance.

If the end user requires to test for a situation-related hygienic-microbiological verification of the dosing device Ecolab suggests the process described in RKI recommendation "requirement on hygiene in cleaning and disinfection of surfaces" (Federal Health Bulletin 2004-47: 51-61, Item 4.5). Kindly follow the following process as described in chapter [8.1 "Microbiological examination of decentralized disinfectant dispensers "](#).

2.3 General Safety Information

- Any connection and repair work on the **DG 1-II** may only be undertaken by authorised and trained experts.
- The **DG 1-II** may only be operated with a supply voltage of 230 V AC.
- Work on electrical equipment must, above all, be undertaken with the mains plug unplugged.
- When changing the disinfectant canister, the warning and safety information on the original disinfectant canister must be observed.
- Safety regulations covering the handling of chemicals must always be observed.

 ATTENTION	<p>At units with Article no. 1024: If using products <u>that contain glucoprotamin</u> (such as INCIDIN PLUS) the O-rings of the metering pump piston and piston rods and the O-rings of the flow monitor must be converted for EPDM. (conversion kit EPDM on request).</p>
	<p>At units with Article no. 102410: If using products <u>that do not contain glucoprotamin</u>, we recommend that you convert the O-rings of the metering pump piston and piston rod and the O-rings of the flow monitor for silicone (conversion kit Silicone on request).</p>
	<p>When using the conversion kit, the provided sticker must be attached to the metering unit in a visible position.</p>

2.4 Maintenance and Repair

Repair and maintenance work must be carried out only by authorized, trained service technicians.

 CAUTION	<p>For maintenance and repair work suitable protective equipment (PSA) has to be worn. Opening covers or removing parts, even if possible without tools, may expose hot wires. Connecting points may also carry an electrical current.</p>
 IMPORTANT	<p>Only original spare and replacement parts may be used for repairs.</p>

2.5 Safety measures (to be performed by the owner/operator)

Attention is drawn to the fact that the owner/operator has the following responsibilities towards his operating and maintenance personnel:

- To provide briefing on the **protective devices** on the equipment;
- Monitor observance of the safety measures;

This operating manual must be retained for future use. The frequency of inspections and control measures must be observed.

The tasks described in this Operating Manual are explained in such a way that:

- any tasks outlined in the chapters Operation and Operating Modes are carried out by a **trained person**,
- the tasks explained in the chapters Transportation, Setup and Assembly, Errors / Causes / Remedies are understood by a **qualified person**.

The chapters Transportation, Setup and Assembly, Errors / Causes / Remedies are only intended for **qualified personnel**. Tasks that are described in these chapters may only be performed by qualified personnel.

2.5.1 Trained person

A person, who has been trained in the tasks assigned to him/her by a skilled person in the possible dangers that improper handling can cause and has been taught about the necessary protective equipment and safety measures.

2.5.2 Qualified person

A person with suitable training, and having suitable education and experience, who is in the position of being able to recognize the risks and knows how to avoid hazards.

The definition is based on EN 60204-1:2007.

2.6 Responsibilities of the Owner

The owner must observe that:

- Personnel safety (Ordinance on Industrial Safety and Health and Accident Prevention Guidelines, Workplace Directives), e.g. operating manuals, also pursuant to §20 GefStoffV, personal protective equipment (PPE)), preventive examinations;
- the safety of the work equipment (protective equipment, instructions for work, process risks and maintenance);
- the product condition (Safety Data Sheets, Register of Hazardous Substances);
- product disposal (the Waste Act);
- material disposal (de-commissioning, Waste Act);
- cleaning (cleaning agents and disposal),
- as well as the latest environmental protection requirements.

In addition, the owner must:

- provide the *personal protective equipment* (PPE): Safety footwear, safety gloves, safety eyewear, protective mask, protective clothing.
- guarantee the training of the operating personnel;
- prepare *work procedures*;

Connections: Before operating the equipment, the owner/operator must ensure that local regulations (e.g. for the electrical connection) were observed during the assembly and commissioning, if this work was carried out by the owner/operator himself.

2.7 System Separation

The integrated system separation according to DIN EN 1717 prevents the intruding of disinfection fluids into the drinking water pipeline in the case of a vacuum in the water supply line.

2.8 Emphases / Numberings

The emphases indicated here have the following meanings:

	CAUTION	This heading is used where imprecise or non-adherence to the operating instructions, prescribed work processes or similar could lead to injuries or accidents.
	ATTENTION	This heading is used where imprecise or non-adherence to the operating instructions, prescribed work processes or similar could lead to damage to the equipment.
	IMPORTANT	Is used when special attention is required when for handling the equipment.
	NOTE	This heading is used to indicate an aspect of particular importance.

✂ Numberings marked with this sign (✂) describe an activity which must only be performed by the installer / user.


2.9 Safety symbols

	DANGER	Dangers due to electrical power are identified with the symbol shown here.
	DANGER	Always wash hands before breaks or at the end of work. Observe the precautionary measures that are usual when handling chemicals. Observe the respective Safety Data Sheet.
	WARNING	When working in areas that are identified by the symbol shown here, it is necessary to wear safety eyewear.
	WARNING	When working in areas that are identified by the symbol shown here, it is necessary to wear safety face protection shield.
	WARNING	When working in areas that are identified by the symbol shown here, it is necessary to wear protective clothing.
	WARNING	When working in areas that are identified by the symbol shown here, it is necessary to wear protective gloves.
	WARNING	When working in areas that are identified by the symbol shown here, it is necessary to wear safety footwear.
	NOTE	Environment symbol identifies environmental safety measures.

2.10 Machine identification

	NOTE	<p><u>The correct information is important for all queries:</u></p> <ul style="list-style-type: none"> • the name • the type • the consecutive serial number of the unit <p>All of this information is required for specific and fast processing. For the English version of DG1-II a corresponding front label can be obtained on request!</p>
---	-------------	--

3 Design

Pos.	Description
1	Cover
2	Mains cable with plug
3	Operation indicator (green LED)
4	Product shortage and fault indicator (red LED)
5	ON/OFF switch (I/O)
6	Bracket
7	Swivel drain outlet
8	Screw (2 x)

Fig. 3.1


Fig. 3.2

Pos.	Description
1	Cover
2	Mains cable with plug
3	Operation indicator (green LED)
4	Product shortage and fault indicator (red LED)
5	ON/OFF switch (I/O)
6	Bracket
7	Swivel drain outlet
8	Screw (2 x)
9	Rear housing cover
10	Coupling nut
11	Control unit
12	Suction pipe
13	Hydraulic motor
14	Eccentric (concentration setting)
15	Mixing chamber
16	Sealing plug
17	Fixing screw
18	Flow control device
19	Dowel Pin
20	Solenoid valve
21	Metering pump
22	Flow interrupter

4 Functional Description and Operation

The **DG1-II** is a mixing equipment for disinfectants. It is used for the production of a ready-to-use disinfectant solution.

The **DG1-II** metering unit includes a double-action piston metering pump driven by a hydraulic motor. The disinfectant is metered proportionally to the water flow rate. By means of an eccentric located at the hydraulic motor, you can set the concentration of the disinfecting solution to a value between 0.25 and 3 %.

IMPORTANT To make sure that the concentration is maintained at its correct value, the minimum tapping quantity must be 1 liter.

4.1 Description of Display Screen and Operating Controls

On the front of the **DG 1-II** is a control panel. It has a membrane keypad and an error indicator (red LED). You use the membrane keypad to operate the unit.

To start or stop the tapping process, press the ON/OFF (I/O) switch. In case of water shortage or a failure in the unit (metering pump defective, flow control device defective, etc.), the metering process is stopped immediately. After eliminating the cause of the trouble, press the (I/O) switch to reset the unit.

Indicator: **Green LED:** Operation/standby indicator (Fig. 3.1, pos. 3)

- Permanent light indicates standby mode
- Flashing light indicates metering mode

Red LED: Product shortage and fault indicator (Fig. 3.1, pos. 4)

- Permanent light indicates a fault
- Flashing light indicates product shortage

The DG 1-II is equipped with the following safety devices:

Water shortage monitoring: The electronic system indirectly monitors the water flow pressure through intelligent evaluation of the pulses emitted by the flow control device.

Flow rate limitation: The maximum quantity of water supplied is controlled by a quantity regulator located inside the solenoid valve.

Product shortage monitoring: The flow control device with its integrated pressure control valve which is fitted upstream from the metering point monitors the flow of the disinfectant product. If a failure of the metering pump, suction line, drive, or level suction pipe causes a product shortage, or if an unsuitable disinfection product is used, the metering unit changes to fault mode. The red fault LED lights up (permanent light).

Suction pipe with integrated level electrodes: If the supply container is empty, the metering unit stops operating and the red fault LED flashes.

Backflow preventer: Upstream from the disinfection metering point, a backflow preventer according to DIN EN 1717 is fitted. This is a safeguard against disinfectant penetrating into the water supply network in case of a backflow.

5 Connection and Installation

The **DG1-II** is designed for the connection of the product reservoir on the right-hand side. If required, the product reservoir can be placed on the left-hand side of the housing. This involves moving the suction lance and the outlet for the power cable to the left-hand side of the housing.

	<p>NOTE</p> <p>The following connection and installation instructions reflect the recommended method. The connection and installation method to be used, however, depends on the actual environmental conditions. All local regulations must be adhered to at all times.</p> <p><u>Housing measurements and connection diagram:</u> see chapter 9.2 „Connection Diagram „Cold Water““ eg. chapter 9.3 „ Connection Diagram „Hot Water““.</p>
---	---

The **DG1-II** is designed to be mounted on a wall, preferably above or close to a sink.

The water supply temperature must not exceed 30°C. In case of higher water temperatures, a thermal mixer valve must be fitted upstream.

The water flow pressure should be at least 0.16 MPa (1.6 bar) and must not exceed 0.7 MPa (7 bar).

If the water flow pressure is > 0.4 MPa (4 bar), we urgently recommend to fit a pressure governor in the water supply line!

The supply voltage is 230 V AC

	<p>NOTE</p> <p>Ex works, the unit is shipped with the suction pipe and mains cable fitted on the right-hand side. If required, the suction pipe and the mains cable can be shifted to the left-hand side of the unit (see chapter 5.4).</p> <p>This operation should be carried out before the unit is mounted to the wall.</p>
---	---

Delivery:

- 1 pce. swivel drain outlet (length: 250 mm)
- 1 pce. angled valve
- 1 set of attachment parts
- 1 pce. drilling jig
- 1 pce. I & O Manual

5.1 Wall Mounting

	<p>ATTENTION Depending on the wall condition, the use of special attachment parts may be necessary.</p>
---	--

	<p>NOTE</p> <p>Be sure to mount the unit at a location high enough to be able to place a drain pan below the swivel drain outlet. The installation height of approx. 130 cm bottom line DG 1, should not overstep (risk of accident when changing the can).</p>
---	--

To mount the DG 1-II proceed as follows:

- ✘ Mark drilling holes for metering unit and bracket using the supplied drilling jig.
- ✘ Drill marked holes..
- ✘ Remove cover (see Fig. 3.1, pos. 1)
for this purpose, unscrew the two cover screws (Fig. 3.1, pos. 8)
Lift the cover to remove it.
- ✘ Mount unit and bracket on the wall.

5.2 Water connection

- ✘ Connect angled valve with its built-in strainer to the wall water supply using a T-union, if necessary.
- ✘ Route the connecting line (suitable type: chromium plated copper pipe, Ø 10 X 1) between the angled valve and the water connector at the bottom of the unit.
- ✘ Attach the supplied swivel drain outlet to the bottom of the unit.


Fig. 5.1

ATTENTION Do not load the swivel spout downwards.

5.3 Electrical Connection

The metering unit may only be connected to a supply voltage of 230 V AC with 50/60 Hz.

For the electrical connection of the **DG 1-II**, all valid national and local regulations must be complied with.

The integrated power transformer of the **DG 1-II** supplies the entire unit with electric power. The unit is protected by a fuse which is located on the PCB (630 mA, slow blow) and by a thermal fuse inside the power transformer.

For the connection to the power supply, the **DG 1-II** is equipped with a mains cable (length: 2.4 m) with an integrated plug. If necessary, the mains cable can be retracted into the housing.

NOTE If required, the mains cable can be shifted to the left-hand side of the housing (see chapter 5.4 „Shifting Suction Pipe and Mains Cable to Left-Hand Side of Housing“).

5.4 Shifting Suction Pipe and Mains Cable to Left-Hand Side of Housing

To shift suction pipe and mains cable to the left-hand side of the housing, proceed as follows:

- ✘ Remove the cover (Fig. 3.1, pos. 1).
For this purpose, unscrew the two cover screws (Fig. 3.1, pos. 8).
Lift the cover to remove it.
- ✘ Detach the hose clamp, unscrew the coupling nut (Fig. 3.2, pos. 10) from the suction pipe, and pull the suction pipe, along with the suction line and the connecting cable, out of the housing.
- ✘ Pull the mains cable bushing out of its holder and retract it into the housing, together with the mains cable.
- ✘ Remove the cover plate at the left-hand side of the housing.
- ✘ Use the cover plate to close the cable opening at the right-hand side of the housing.
- ✘ Pull the mains cable through the left-hand opening of the housing and attach the cable bushing. Make sure that the cable is not damaged and cannot get in touch with any sharp-edged parts (leave the internal mains cable attachments in place).
- ✘ Fit the suction pipe to the left-hand side of the housing

6 Settings and Start-Up

6.1 Start-Up

**NOTE**

All item numbers in this section refer to the exploded view unless otherwise stated.

**ATTENTION**

If products are used which contain glucoprotamine (e.g. INCIDIN PLUS), the O-rings on the piston and on the piston rod need to be replaced by EPDM O-rings.

After connection and installation work proceed as follows:

- ✘ Open the angled valve (G ½) to its full open position, and check it for tightness.
- ✘ Connect the disinfectant supply container as follows:
 - ✘ swivel suction pipe to front position and insert it into opening of supply container,
 - ✘ place the supply container onto the bracket.
- ✘ Bleed the metering pump as follows:
 - ✘ unscrew the flow control device (see Fig. 3.2, pos. 18) (black coupling nut) and place it in a suitable jar,
 - ✘ Close the opening of the mixing chamber (see Fig. 3.2, pos. 16) by inserting the sealing plug (see Fig. 3.2, pos. 15),
 - ✘ unscrew the knurled nut (see Fig. 6.1, pos. 3) and withdraw the eccentric (see Fig. 6.1, pos. 1) from its holder located on the hydraulic motor (see Fig. 6.1, pos. 6),
 - ✘ move the connecting rod including the piston rod (see Fig. 6.1, pos. 7) forth and back until disinfectant runs out from the flow control device and no more bubbles are visible in the hoses,
 - ✘ remove the sealing plug and screw in the flow control device again with the VA pins pointing downwards,
 - ✘ plug in the eccentric back in and tighten the knurled nut.
- ✘ Set the desired concentration as described in chapter [6.2](#) „Setting the Concentration“ and check it as described in chapter [6.2.1](#) „Checking the Concentration Setting“.
- ✘ Check the safety equipment as described in chapter [6.2.2](#) „Checking the Safety Equipment“.
- ✘ Check all lines and components which come in contact with water and/or disinfectant for tightness.

6.2 Setting the Concentration

Two concentration ranges can be set:

- Concentration range 1: 0,25 - 1 % (white eccentric)
- Concentration range 2: 1,5 - 3 % (blue eccentric)

Ex works, the unit is set to a concentration of 0.5 %.


NOTE

If you plan to use concentration range 2, the white eccentric must be replaced by the supplied blue eccentric. For this purpose, remove the knurled nut (see Fig. 6.1, pos. 3) and the graduated dial (see Fig. 6.1, pos. 4).

ATTENTION When setting the concentration, allow for a safety margin of + 10 %.

Pos.	Description
1	Eccentric
2	Red mark
3	Knurled nut
4	Graduated dial
5	Metering pump
6	Drive unit (hydraulic motor)
7	Piston rod

Fig. 6.1

To set the nominal concentration, proceed as follows:

- ✘ Unscrew knurled nut (see pos. 3) from eccentric (see pos. 1) by approx. 5 mm.
- ✘ Pull graduated dial (see pos. 4) from indented part of eccentric and hold it in place.
- ✘ Preset the concentration: for this purpose, turn the eccentric so that the red mark aligns with the desired concentration value on the graduated dial.
- ✘ Insert the graduated dial back into the indented part of the eccentric and tighten the knurled nut.
- ✘ Check the concentration as described in chapter 6.2.1 „Checking the Concentration Setting“. If it deviates from the nominal value, repeat the setting procedure.
- ✘ Align the arrow of the red mark on the eccentric with the appropriate concentration value on the graduated dial. For this purpose, loosen the locking screw on the eccentric.

6.2.1 Checking the Concentration Setting

Required auxiliaries: Bucket with a minimum capacity of 10 l to collect the water

Measuring jars to collect the disinfectant:

up to 100 ml for concentrations of 1 % max.

up to 300 ml for concentrations of 1 % to 3 %

In the following example which shows how to determine the nominal concentration, the quantity settings given in millimeters refer to a tapped water quantity of exactly 8 liters.

Formulas for setting the concentration:

Explanantion:

K = Concentration (%)
 VD = Volume of disinfectant (ml)
 VW = Volume of water (ml)

Formula:

$$K = \frac{VD \times 100}{VW + VD} \quad VD = \frac{K \times VW}{100 - K}$$

Example: Desired concentration = 0.5 %, plus 10 % buffer = 0.55 %

Disinfectant quantity to be set according to the above formula:

$$VD = \frac{K \times VW}{100 - K} = \frac{0,55 \times 8.000}{100 - 0,55} = \mathbf{44,24 \text{ ml}}$$

To check the concentration setting, proceed as follows:

- ✘ Unscrew the flow control unit from the mixing chamber and seal the mixing chamber opening with the supplied sealing plug.
- ✘ Insert the flow control unit into the provided measuring jar and place a bucket below the swivel drain outlet.
- ✘ Press the ON/OFF (I/O) switch and tap approx. 2-3 liters of water.
- ✘ Empty the bucket and the measuring jar.
- ✘ Insert the flow control unit back into the measuring jar and place the bucket below the swivel drain outlet again.
- ✘ Start the unit by pressing the ON/OFF (I/O) switch.
- ✘ Tap exactly 8 liters of water. Switch off the unit (press the ON/OFF (I/O) switch).
- ✘ Read the disinfectant quantity collected in the measuring jar.
- ✘ Calculate the concentration according to the above formula. If it deviates from the desired concentration (plus buffer), adjust the setting as described in chapter [6.2](#) „Setting the Concentration“.
- ✘ Remove sealing plug from mixing chamber opening and screw flow control device in again.

6.2.2 Checking the Safety Equipment

A check of the safety equipment obligatory includes the following test procedures:

- functional test of the level suction pipe
- functional test of the flow control device

6.2.3 Functional Check of the Level Suction Pipe

For the functional test of the suction pipe "empty" signal device, proceed as follows:

- ✘ Pull suction pipe out of product supply container.
- ✘ After approx. 3 seconds, unit should change to fault mode (red fault indicator flashes).

6.2.4 Functional Check of the Flow Control Device

For the functional check proceed as follows:

- ✘ Remove the cover.
- ✘ Unscrew the knurled nut (see Fig. 6.1, pos. 3) from the eccentric and remove the complete eccentric from its holder on the hydraulic motor.
- ✘ Press the ON/OFF (I/O) switch. Approx. 10 seconds after the unit has been started, it should change to fault mode (the red LED is permanently ON).
- ✘ Refit the eccentric.
- ✘ Press the ON/OFF (I/O) switch to reset the unit.
- ✘ Press the ON/OFF (I/O) switch again (unit starts metering).
- ✘ Slowly reduce the supply of fresh water to reach 0 MPa (0 bar) of pressure. When the pressure falls below 1.4 MPa (1.4 bar) the unit changes to fault mode (red LED is permanently ON).
- ✘ Reopen the fresh water supply.
- ✘ Press the ON/OFF (I/O) switch to reset the unit.
- ✘ Reposition the cover.

6.2.5 Changing over to a New Product

ATTENTION

When changing over to a new product, make sure to thoroughly flush the suction pipe, the metering pump, and all lines in contact with the product before you connect the new supply container.

If products are used which contain glucoprotamine (e.g. INCIDIN PLUS), the O-rings on the piston and on the piston rod need to be replaced by EPDM O-rings.

If products are used which don't contain glucoprotamine, the O-rings on the piston and on the piston rod need to be replaced by SILICONE O-rings (see chapter 2.3 „General Safety Information“).

To change over to a new product, proceed as follows:

- ✘ Fill an appropriate vessel with water and immerse the suction pipe
- ✘ Flush suction pipe, metering pump, and the lines in contact with the product as follows:
 - unscrew the flow control device (see Fig. 3.2, pos 12), insert the sealing plug, and place the flow control device in an appropriate vessel,
 - unscrew the knurled nut (see Fig. 6.1, pos 3) from the eccentric and pull the eccentric from its holder on the hydraulic motor,
 - manually move the piston rod (see Fig. 6.1, pos 7) forth and back until clear water without product residues is running from the flow control device.
- ✘ Insert the suction pipe into the new supply container.
- ✘ Bleed metering pump until product running out of flow control device is bubble-free (see chapter 6.1 „Start-Up“)
- ✘ Push the eccentric back onto its holder and tighten the knurled nut.
- ✘ Remove the sealing plug and screw the flow control device in again.
- ✘ Check the unit for leaks and malfunctions.
- ✘ Set the concentration as described in chapter 6.2 „Setting the Concentration“ or check it as described in chapter 6.2.1 „Checking the Concentration Setting“.

7 Troubleshooting

 **NOTE** Each time you stop the unit and/or take corrective action to eliminate a fault, make sure to reset the unit by pressing the ON/OFF (I/O) switch.

Failure	Cause	Corrective Action
Control unit does not work (no LED lights up, pressing a button has no effect)	No supply voltage available	Restore power supply
	Fuse of control unit defective	Replace fuse
	Control PCB defective	Replace control PCB
	Thermal fuse of power transformer blown	Let metering unit cool down
	Power transformer defective	Replace power transformer
Pressing the ON/OFF (I/O) switch does not start the unit; fault indicator <u>is flashing</u>	Disinfectant container empty	Connect new supply container
	Suction pipe soiled, suction pipe defective	Clean bottom valve of suction pipe
	Connecting cable of suction pipe defective	Replace suction pipe
Fault indicator <u>permanently ON</u>	Metering pump leaks or defective	Replace metering pump
	Suction and metering line: Hose clamps leak; hose defective; angled valve closed	Tighten hose clamps, replace hose, if necessary open angled valve
	Water flow pressure < 0.16 MPa bar (1.6 bar)	Ensure that water flow pressure is > 0.16 MPa (1.6 bar)
	Strainer clogged	Clean strainer, replace it, if necessary
	Flow limiter encrusted with lime or clogged	Clean flow limiter, replace it, if necessary
	Hydraulic motor leaks	Replace axial face seal
	Hydraulic motor blocked	Replace hydraulic motor
	Coil of solenoid valve defective	Replace coil
Unit switched OFF; swivel drain outlet leaks	Solenoid valve soiled or defective	Clean solenoid valve, replace defective parts
Concentration does not reach the set value	Metering pump: Suction and pressure valves soiled or defective	Clean valves, replace, if necessary
	Suction pipe: Bottom valve soiled or defective	Clean bottom valve, replace, if necessary
	Water feed pipe soiled	Clean strainers in valve and in fresh water supply pipe, replace, if necessary
	Hydraulic motor defective	Replace hydraulic motor
Product shortage monitoring device does not switch the unit off when the product container is empty	Suction pipe electrodes soiled or defective, connecting cable damaged	Replace suction pipe
Unit cannot be reset	Control PCB defective	Replace control PCB

8 Maintenance and Repair Information

CAUTION

Suitable protective equipment must be worn when servicing or repairing parts which come into contact with disinfectant. The safety regulations which apply for handling chemicals must always be observed.

ATTENTION

To ensure operational safety and correct operation of the **DG 1-II**, the specified maintenance intervals must be respected. The interval between the installation and the first maintenance session or the period between two consecutive maintenance sessions must not exceed 12 months.

8.1 Microbiological examination of decentralized disinfectant dispensers

- ✘ Unscrew nut from S- spout of unit and remove aerator.
- ✘ The outside of outlet has to be to disinfected with an alcoholic surface disinfectant (Incidin Liquid, Incides N) or to treated with a flame
- ✘ After flow of approx. 1 liter collect at least minimum 1 liter of disinfectant solution from the device in a sterile sampling container (Note: 1L is the minimum amount to be drawn, see also note on the device).
- ✘ Neutralizing agents must only be added after recommended contact time for the relevant use concentration of disinfectant against bacteria and yeast (see product label).
- ✘ Appropriate neutralizers (chemical composition and concentration) can be selected and validated according to EN 13727 ort the Standard Methods of DGHM/VAH or be asked from the disinfectant manufacturer
- ✘ Neutralizing agent must be sterile
- ✘ Sample has to be investigated in a microbiological laboratory by filtration method for total microbial count. Ideally, this value is determined per 1 ml, per 10ml and per 100 ml
- ✘ A refrigerated transport (4 ° C) in lightproof containers is desirable; it is mandatory to transport the samples within 3 hours of sample collection.
- ✘ The treatment of the sample should be completed 2 to 3 hours after taking out.
- ✘ During the assessment process of the results, identification of the micro-organism is also recommended and not only the microbial count to allow assessment of results with regard to the efficacy spectrum of the disinfectant dosed.

8.2 General maintenance work

After completion of the maintenance work, perform the following tasks:

- ✘ **Care of internal parts:** clean bottom plate of unit and remove any lime residues.
- ✘ **Care of external parts:** Clean the cover with a damp cloth, rub it dry, and then treat it with product for stainless steel care (e.g. Chromol).
- ✘ **Bracket::** Clean surface where you place supply container

To perform maintenance work on the **DG 1-II**, proceed as specified in the following Maintenance Table (see chapter [8.2.1](#) „Maintenance Table“).

8.2.1 Maintenance Table

Component/Maintenance Task - Check Criterion-	Work Instructions/References
<u>Flow control device:</u> – Functional test – Tightness of pressure control valve	See chapter 6.2.4 „Functional Check of the Flow Control Device“ If disinfectant is running from the flow control device, unscrew it from the threaded union ⇒ replace flow control device
<u>Checking the concentration setting:</u>	See chapter 6.2.1 „Checking the Concentration Setting“
<u>Eccentric:</u> Functional test	To check the eccentric for correct operation, set the concentration to three different values and gauge them by liters. In case of deviation, replace eccentric, connecting rod, and – if necessary – complete metering pump
<u>Metering pump:</u> Leak test	With the unit running, watch for concentrate escaping at the piston rod. If yes ⇒ replace piston and O-rings or replace the complete metering pump, if necessary
<u>Hydraulic motor:</u> Performance test	With the unit switched on, grasp the eccentric. If it is possible to block the hydraulic motor without applying force: Hydraulic motor defective ⇒ replace hydraulic motor
<u>Backflow preventer:</u> Leak test of DB-type flow interrupter	Replace rubber diaphragm or flow interrupter, if necessary
<u>Suction pipe:</u> Functional test	See chapter 6.2.3 „Functional Check of the Level Suction Pipe“
<u>Suction pipe:</u> Leak test	Immerse suction pipe end and electrodes in a beaker filled with disinfectant; if the metering pump sucks air, the suction pipe hose is defective ⇒ replace hose or suction pipe, if necessary
<u>Strainer: (Fig. 3.3, pos. 7)</u> One or more strainers soiled Solenoid valve strainer	Remove insets from the affected strainers and clean them. Remove the coupling nut from the slanted-seat filter.

8.3 Recommendations for the replacement of wearing parts

Component	Work step	Maintenance interval
Solenoid valve:	Exchange of the filter	annual
Connectinghose/Water -separation & outlet		annual
System separation (pipe interrupter DB)	Exchange of the diaphragm	every 3 years

8.4 Repair Instructions

Repair instructions including spare parts lists and exploded views you can get on request.

	<p>CAUTION</p> <p>For any work on parts which come in contact with the product, the corresponding hazard information and safety instructions must be observed at all times. Always wear suitable protective equipment.</p> <p>Before carrying out any repair work, make sure to disconnect the metering unit from the mains. To do so, pull the mains plug out of the wall outlet.</p> <p>Before performing any repair work on parts which come in contact with the product and/or with water (e.g. solenoid valve), make sure that the angled valve is closed and that the water pressure inside the unit is down. To do so, press the ON/OFF (I/O) switch twice.</p>
	<p>ATTENTION</p> <p>Only original spare parts may be used! (see chapter 1.3 „Scope of guarantee“).</p>

9 Technical Data

9.1 General Data


Supply voltage:	AC 230 V
Power consumption:	15 W
Protection of metering unit/control unit:	IP 24 / IP 65
Protection class:	II according to EN 60335-1
Superposed fuse:	max. 16 A
Operating voltage:	AC 24 V
Internal fuse protection:	Microfuse 5 x 20, 630 mA, slow blow
Water connection:	Conical coupling NW 10 (for copper pipe, D. 10x1)
Water supply temperature:	max. 30°C
Flow pressure::	min. 0.16 MPa (1.6 bar) - max. 0.7 MPa (7 bar)
Metering capacity:	approx. 850 l/h (dilution for normal operation)
Concentration setting:	White eccentric: 0.25 % - 1 %, Blue eccentric: 1.5 % - 3 %
Metering tolerance:	Max. + 6 %
Backflow preventer:	DIN EN 1717 (Typ DB)
Mounting dimensions:	Width = 500 mm (with bracket), High =425 mm, Depth=200 mm
Weight:	13,5 kg

**NOTE**

Due to the company's policy of continuous improvements to its products, the specifications may be changed without notice!

9.2 Connection Diagram

Fig. 9.1

Pos.	Connections
A	Water for Cu pipe Ø 10 x 1
B	Power cable AC 230 V
C	Swivel spout

Swivel spout	L	H
202526	150	138
202513	200	149
Standard 202514	250	151
202515	300	153
202416	280-480	145
202415	350-600	145

9.3 Connection Material and Accessories

Connection Material - Designation
T-Union G 1/2, MS chromium-plated brass
T-Union R 3/4 / R 1/2 / R 3/4, MS chromium-plated brass
Copper pipe 10 x 1 chromium-plated brass
Angled valve G 1/2 a x D 10 MS chromium-plated brass [supplied]
Accessories - Designation
Adapter 3/4 " internal thread R1/2 external thread, chromium-plated brass
Swivel drain outlet, length 150 mm
Swivel drain outlet, length 200 mm
Swivel drain outlet, length 300 mm
Thermal mixer valve unit including screw coupling
Conversion kit for restoring the unit to a clearance pipe segment configuration (only for to low water pressure up to 0,1 MPa (1 bar)
Bracket for 10L-Can

NOTE If necessary please ask for part numbers together with the spare part drawing separately!

10 Dismantling / Disposal

Dismantling

The **DG1-II** may only be dismantled by specialist personnel. Before commencing dismantling, ensure that the **DG1-II** has been fully isolated from the power supply. The device must be flushed through with care in order to get rid of any residual chemicals.

Disposal


The **DG1-II** is predominantly made of steel (and also of aluminium to a certain extent) (apart from the electrical equipment) and is to be disposed of in accordance with the local environmental regulations prevailing **at that time**.

Depending on the condition, existing regulations and with due regard for current requirements, dispose of it, for example, as

electrical waste (printed circuit boards), plastics (housing), metal plates, steel, copper, aluminium (separate after sorting). Prior to disposal, all parts which are in contact with media must be decontaminated.

Oils, solvents, detergents and contaminated cleaning tools (brushes, cloths, etc.) must be disposed of in compliance with local requirements, in accordance with the prevailing waste code and with due attention to the notes contained in the manufacturers' safety data sheets.

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Appendix 1 to WI-EE-RD&E-04

03 / 2012

	<p>NOTE</p> <p>This declaration of conformity corresponds to the status of this operating instruction. Should in the meantime one of those standards be subject of change, you can download the latest version here under: http://www.ecolab-engineering.de/uploads/media/CE_1024ff_DG1.pdf</p>	
---	--	---

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© Copyright **Ecolab Engineering GmbH**, 2016
Alle Rechte vorbehalten. *All rights reserved*
ECOLAB Engineering GmbH Postfach 11 64 D-83309 Siegsdorf
der Firma **Ecolab Engineering GmbH** gestattet.